

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502631

UBXN2B Mouse Monoclonal Antibody [Clone ID: OTI1F7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1F7

Applications: FC, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human UBXN2B (NP_001071087) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.41 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 36.9 kDa

Gene Name: UBX domain protein 2B


Database Link: NP 001071087 Entrez Gene 137886 Human

Synonyms: p37






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY UBXN2B ([RC219137], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-UBXN2B. Positive lysates [LY421462] (100ug) and [LC421462] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-UBXN2B mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502631)







Immunohistochemical staining of paraffinembedded Human lymphoma tissue using anti-UBXN2B mouse monoclonal antibody. (Heatinduced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502631)



HEK293T cells transfected with either [RC219137] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-UBXN2B antibody (TA502631), and then analyzed by flow cytometry.



Flow cytometric Analysis of Hela cells, using anti-UBXN2B antibody (TA502631), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-UBXN2B antibody (TA502631), (Red), compared to a nonspecific negative control antibody ([TA50011]), (Blue).