

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502755

RBP1 Mouse Monoclonal Antibody [Clone ID: OTI2H3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2H3

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human RBP1 (NP_002890) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1.2 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 15.7 kDa

Gene Name: retinol binding protein 1

Database Link: NP 002890 Entrez Gene 5947 Human

Background: This gene encodes the carrier protein involved in the transport of retinol (vitamin A alcohol)


from the liver storage site to peripheral tissue. Vitamin A is a fat-soluble vitamin necessary for growth, reproduction, differentiation of epithelial tissues, and vision. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq]

Synonyms: CRABP-I; CRBP; CRBP1; CRBPI; RBPC






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RBP1 ([RC229119], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RBP1. Positive lysates [LY401017] (100ug) and [LC401017] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (10ug) from 9 Human tissue by using anti-RBP1 monoclonal antibody at 1:200 (1: Testis; 2: Omentum; 3: Uterus; 4: Breast; 5: Brain; 6: Liver; 7: Ovary; 8: Thyroid gland; 9: colon).






Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-RBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502755)




Immunohistochemical staining of paraffinembedded Human Ovary tissue within the normal limits using anti-RBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502755)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-RBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502755)






Immunohistochemical staining of paraffinembedded Human endometrium tissue within the normal limits using anti-RBP1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502755)




Anti-RBP1 mouse monoclonal antibody (TA502755) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY RBP1 ([RC214515]).




Immunofluorescent staining of HeLa cells using anti-RBP1 mouse monoclonal antibody (TA502755).






HEK293T cells transfected with either [RC214515] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-RBP1 antibody (TA502755), and then analyzed by flow cytometry.



Flow cytometric Analysis of Hela cells, using anti-RBP1 antibody (TA502755), (Red), compared to a nonspecific negative control antibody, (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-RBP1 antibody (TA502755), (Red), compared to a nonspecific negative control antibody, (Blue).