

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA502891

eRF1 (ETF1) Mouse Monoclonal Antibody [Clone ID: OTI2E3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2E3
Applications: FC, WB

Recommend Dilution: WB 1:500, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ETF1 (NP_004721) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.97 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 48.9 kDa

Gene Name: eukaryotic translation termination factor 1

Database Link: NP 004721 Entrez Gene 2107 Human

Background: This gene encodes a class-1 polypeptide chain release factor. The encoded protein plays an

essential role in directing termination of mRNA translation from the termination codons UAA,

UAG and UGA. This protein is a component of the SURF complex which promotes

degradation of prematurely terminated mRNAs via the mechanism of nonsense-mediated mRNA decay (NMD). Alternate splicing results in coding and non-coding transcript variants. Pseudogenes of this gene are found on chromosomes 5, 7 and X. [provided by RefSeq, Jan


2012]

Synonyms: D5S1995; ERF; ERF1; RF1; SUP45L1; TB3-1







Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ETF1 ([RC220387], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ETF1. Positive lysates [LY401492] (100ug) and [LC401492] (20ug) can be purchased separately from OriGene.



Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-ETF1 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human) (1:200).



Flow cytometric Analysis of Hela cells, using anti-ETF1 antibody (TA502891), (Red), compared to a nonspecific negative control antibody, (Blue).



Flow cytometric Analysis of Jurkat cells, using anti-ETF1 antibody (TA502891), (Red), compared to a nonspecific negative control antibody, (Blue).