

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA502935**

## CDK2 Mouse Monoclonal Antibody [Clone ID: OTI2D9]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI2D9

Applications: FC, IF, WB

**Recommend Dilution:** WB 1:500~2000, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human, Monkey, Rat, Dog

Host: Mouse IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human CDK2 (NP\_001789) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.78 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 33.7 kDa

**Gene Name:** cyclin dependent kinase 2

Database Link: NP 001789 Entrez Gene 362817 RatEntrez Gene 100855704 DogEntrez Gene 711002

MonkeyEntrez Gene 1017 Human

**Background:** The protein encoded by this gene is a member of the Ser/Thr protein kinase family. This

protein kinase is highly similar to the gene products of S. cerevisiae cdc28, and S. pombe cdc2. It is a catalytic subunit of the cyclin-dependent protein kinase complex, whose activity is restricted to the G1-S phase, and essential for cell cycle G1/S phase transition. This protein associates with and regulated by the regulatory subunits of the complex including cyclin A or E, CDK inhibitor p21Cip1 (CDKN1A) and p27Kip1 (CDKN1B). Its activity is also regulated by its protein phosphorylation. Two alternatively spliced variants and multiple transcription

initiation sites of this gene have been reported. [provided by RefSeq]

**Synonyms:** CDKN2; p33(CDK2)

**Protein Families:** Druggable Genome, Protein Kinase







**Protein Pathways:** 

Cell cycle, Oocyte meiosis, p53 signaling pathway, Pathways in cancer, Progesterone-mediated oocyte maturation, Prostate cancer, Small cell lung cancer


## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY CDK2 ([RC200494], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-CDK2. Positive lysates [LY419741] (100ug) and [LC419741] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-CDK2 monoclonal antibody.




Western blot analysis of extracts (35ug) from 14 cell lines by using anti-CDK2 monoclonal antibody (1:500).






Equivalent amounts of cell lysates (10 ug per lane) of wild-type 293T cells (WT, Cat# LC810293T) and CDK2-Knockout 293T cells (KO, Cat# [LC810794]) were separated by SDS-PAGE and immunoblotted with anti-CDK2 monoclonal antibody TA502935, (1:2000). Then the blotted membrane was stripped and reprobed with anti-HSP90AB1 antibody ([TA500494]) as a loading control.




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-CDK2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA502935)




Anti-CDK2 mouse monoclonal antibody (TA502935) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY CDK2 ([RC200494]) (1:100).






Immunofluorescent staining of Hela cells using anti-CDK2 mouse monoclonal antibody (TA502935, green, upper left; merged, upper right) or Isotype control (merged, lower right). Cell nuclei were stained with DAPI (blue, lower left) (1:100).




Immunofluorescent staining of 293T cells using anti-CDK2 mouse monoclonal antibody (TA502935, green, upper left; merged, upper right) or Isotype control (merged, lower right). Cell nuclei were stained with DAPI (blue, lower left) (1:100).




HEK293T cells transfected with either [RC200494] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-CDK2 antibody (TA502935), and then analyzed by flow cytometry (1:100).





Flow cytometric Analysis of Hela cells, using anti-CDK2 antibody (TA502935), (Red), compared to a nonspecific negative control antibody, (Blue) (1:100).



Flow cytometric Analysis of Jurkat cells, using anti-CDK2 antibody (TA502935), (Red), compared to a nonspecific negative control antibody, (Blue) (1:100).