

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503185

CRALBP (RLBP1) Mouse Monoclonal Antibody [Clone ID: OTI1H5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1H5

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human

Host: Mouse Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human RLBP1(NP_000317) produced in HEK293

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.76 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 36.3 kDa

Gene Name: retinaldehyde binding protein 1

Database Link: NP 000317 Entrez Gene 6017 Human

Background: The protein encoded by this gene is a 36-kD water-soluble protein which carries 11-cis-

retinaldehyde or 11-cis-retinal as physiologic ligands. It may be a functional component of the visual cycle. Mutations of this gene have been associated with severe rod-cone dystrophy, Bothnia dystrophy (nonsyndromic autosomal recessive retinitis pigmentosa) and retinitis

punctata albescens. [provided by RefSeq]


Synonyms: CRALBP

Protein Families: Druggable Genome





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RLBP1 ([RC201136], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RLBP1. Positive lysates [LY424794] (100ug) and [LC424794] (20ug) can be purchased separately from OriGene.