

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503303

DOK2 Mouse Monoclonal Antibody [Clone ID: OTI1G9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1G9

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DOK2(NP_003965) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.36 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 45.2 kDa

Gene Name: docking protein 2

Database Link: NP 003965 Entrez Gene 9046 Human

Background: The protein encoded by this gene is constitutively tyrosine phosphorylated in hematopoietic

progenitors isolated from chronic myelogenous leukemia (CML) patients in the chronic phase.

It may be a critical substrate for p210(bcr/abl), a chimeric protein whose presence is

associated with CML. This encoded protein binds p120 (RasGAP) from CML cells. [provided by


RefSeq1

Synonyms: p56DOK; p56dok-2 **Protein Families:** Druggable Genome





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DOK2 ([RC207621], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DOK2. Positive lysates [LY418321] (100ug) and [LC418321] (20ug) can be purchased separately from OriGene.