

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA503452

Asparagine synthetase (ASNS) Mouse Monoclonal Antibody [Clone ID: OTI1A10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1A10
Applications: IF, WB

Recommend Dilution: WB 1:2000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ASNS(NP_597680) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 64.2 kDa

Gene Name: asparagine synthetase (glutamine-hydrolyzing)

Database Link: NP 597680 Entrez Gene 440 Human

Background: The protein encoded by this gene is involved in the synthesis of asparagine. This gene

complements a mutation in the temperature-sensitive hamster mutant ts11, which blocks

progression through the G1 phase of the cell cycle at nonpermissive temperature. Alternatively spliced transcript variants have been described for this gene. [provided by

RefSeq1

Synonyms: ASNSD; TS11


Protein Families: Druggable Genome

Protein Pathways: Alanine, aspartate and glutamate metabolism, Metabolic pathways, Nitrogen metabolism






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ASNS ([RC215380], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ASNS. Positive lysates [LY403336] (100ug) and [LC403336] (20ug) can be purchased separately from OriGene.



Anti-ASNS mouse monoclonal antibody (TA503452) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY ASNS ([RC215380]).