

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA504245

C1S Mouse Monoclonal Antibody [Clone ID: OTI5F2]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5F2

Applications: FC, IF, IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human C1S(NP_001725) produced in HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 0.71 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 74.8 kDa

Gene Name: complement C1s

Database Link: NP 001725 Entrez Gene 716 Human

Background: This gene encodes a serine protease, which is a major constituent of the human complement

subcomponent C1. C1s associates with two other complement components C1r and C1q in order to yield the first component of the serum complement system. Defects in this gene are

the cause of selective C1s deficiency. [provided by RefSeq]

Synonyms: FLJ44757


Protein Families: Druggable Genome, Protease

Protein Pathways: Complement and coagulation cascades, Systemic lupus erythematosus







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY C1S ([RC218963], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-C1S. Positive lysates [LY400655] (100ug) and [LC400655] (20ug) can be purchased separately from OriGene.




Anti-C1S mouse monoclonal antibody (TA504245) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY C1S ([RC218963]).




Flow cytometric Analysis of Jurkat cells, using anti-C1S antibody (TA504245), (Red), compared to a nonspecific negative control antibody, (Blue).






Flow cytometric Analysis of Hela cells, using anti-C1S antibody (TA504245), (Red), compared to a nonspecific negative control antibody, (Blue).




Immunohistochemical staining of paraffinembedded Carcinoma of Human bladder tissue using anti-C1S mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504245)




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-C1S mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504245)






Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-C1S mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504245)





Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-C1S mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504245)



Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-C1S mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504245)







Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-C1S mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504245)

Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-C1S mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA504245)