

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA504404

AFG3L2 Mouse Monoclonal Antibody [Clone ID: OTI2D3]

Product data:

Applications:

Product Type: Primary Antibodies

WB

Clone Name: OTI2D3

Recommend Dilution: WB 1:500

Reactivity: Human

Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-250 of human

AFG3L2(NP_006787) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 88.3 kDa

Gene Name: AFG3 like matrix AAA peptidase subunit 2

Database Link: NP 006787 Entrez Gene 10939 Human

Background: This gene encodes a protein localized in mitochondria and closely related to paraplegin. The

paraplegin gene is responsible for an autosomal recessive form of hereditary spastic paraplegia. This gene is a candidate gene for other hereditary spastic paraplegias or

neurodegenerative disorders. [provided by RefSeq, Jul 2008]


Synonyms: SCA28; SPAX5

Protein Families: Druggable Genome, Protease, Transmembrane





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AFG3L2 ([RC209058], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AFG3L2. Positive lysates [LY416420] (100ug) and [LC416420] (20ug) can be purchased separately from OriGene.