

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA504443**

## MIF4GD Mouse Monoclonal Antibody [Clone ID: OTI4C11]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4C11
Applications: FC, IF, WB

Recommend Dilution: WB 1:2000, IF 1:100, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MIF4GD(NP\_065730) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

**Concentration:** 0.78 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 29.2 kDa

**Gene Name:** MIF4G domain containing

Database Link: NP 065730 Entrez Gene 57409 Human

**Background:** This gene encodes a protein which interacts with the N-terminus of the stem-loop binding

protein (SLBP) and the 3' end of histone mRNA. This interaction facilitates the activation of histone mRNA translation. Alternative splicing results in multiple transcript variants encoding


different isoforms. [provided by RefSeq]

Synonyms: AD023; MIFD; SLIP1







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MIF4GD ([RC207489], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MIF4GD. Positive lysates [LY402799] (100ug) and [LC402799] (20ug) can be purchased separately from OriGene.




Anti-MIF4GD mouse monoclonal antibody (TA504443) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MIF4GD ([RC207489]).



Immunofluorescent staining of HeLa cells using anti-MIF4GD mouse monoclonal antibody (TA504443).





HEK293T cells transfected with either [RC207489] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-MIF4GD antibody (TA504443), and then analyzed by flow cytometry.