

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA505031

ARA9 (AIP) Mouse Monoclonal Antibody [Clone ID: OTI1B6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1B6
Applications: FC, WB

Recommend Dilution: WB 1:2000, FLOW 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human AIP(NP_003968) produced in HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 37.5 kDa

Gene Name: aryl hydrocarbon receptor interacting protein

Database Link: NP 003968 Entrez Gene 9049 Human

Background: The protein encoded by this gene is a receptor for aryl hydrocarbons and a ligand-activated

transcription factor. The encoded protein is found in the cytoplasm as part of a multiprotein complex, but upon binding of ligand is transported to the nucleus. This protein can regulate the expression of many xenobiotic metabolizing enzymes. Also, the encoded protein can bind

specifically to and inhibit the activity of hepatitis B virus. [provided by RefSeq]


Synonyms: ARA9; FKBP16; FKBP37; SMTPHN; XAP-2; XAP2

Protein Families: Druggable Genome, Transcription Factors






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY AIP ([RC211157], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-AIP. Positive lysates [LY401308] (100ug) and [LC401308] (20ug) can be purchased separately from OriGene.



HEK293T cells transfected with either [RC211157] overexpress plasmid (Red) or empty vector control plasmid (Blue) were immunostained by anti-AIP antibody (TA505031), and then analyzed by flow cytometry.