

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA505709**

### ATXN7L1 Mouse Monoclonal Antibody [Clone ID: OTI3G10]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI3G10

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ATXN7L1 (NP\_689962) produced in

HEK293T cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 16 kDa


**Gene Name:** ataxin 7 like 1

Database Link: NP 689962 Entrez Gene 222255 Human

Synonyms: ATXN7L4



# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ATXN7L1 ([RC206079], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ATXN7L1. Positive lysates [LY407300] (100ug) and [LC407300] (20ug) can be purchased separately from OriGene.