

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA505914

RNF39 Mouse Monoclonal Antibody [Clone ID: OTI4D3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4D3
Applications: IF, WB

Recommend Dilution: WB 1:2000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human RNF39(NP_079512) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 45.3 kDa

Gene Name: ring finger protein 39


Database Link: NP 079512 Entrez Gene 80352 Human

Synonyms: HZF; HZFW; LIRF


Protein Families: Druggable Genome



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RNF39 ([RC219332], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RNF39. Positive lysates [LY410747] (100ug) and [LC410747] (20ug) can be purchased separately from OriGene.



Anti-RNF39 mouse monoclonal antibody (TA505914) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY RNF39 ([RC219332]).