

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA506034**

### Junctional Adhesion Molecule 1 (F11R) Mouse Monoclonal Antibody [Clone ID: OTI6E11]

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI6E11
Applications: IF, WB

Recommend Dilution: WB 1:4000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human F11R(NP\_653087) produced in HEK293T

cell

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 29.8 kDa **Gene Name:** F11 receptor

Database Link: NP 653087 Entrez Gene 50848 Human

**Background:** Tight junctions represent one mode of cell-to-cell adhesion in epithelial or endothelial cell

sheets, forming continuous seals around cells and serving as a physical barrier to prevent solutes and water from passing freely through the paracellular space. The protein encoded by this immunoglobulin superfamily gene member is an important regulator of tight junction assembly in epithelia. In addition, the encoded protein can act as (1) a receptor for reovirus, (2) a ligand for the integrin LFA1, involved in leukocyte transmigration, and (3) a platelet receptor. Multiple 5' alternatively spliced variants, encoding the same protein, have been identified but their biological validity has not been established. [provided by RefSeq, Jul 2008]

**Synonyms:** CD321; JAM; JAM-1; JAM-4; JAM1; JAMA; JCAM; KAT; PAM-1

**Protein Families:** Druggable Genome, Transmembrane







**Protein Pathways:** 

Cell adhesion molecules (CAMs), Epithelial cell signaling in Helicobacter pylori infection, Leukocyte transendothelial migration, Tight junction

## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY F11R ([RC200004], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-F11R. Positive lysates [LY403392] (100ug) and [LC403392] (20ug) can be purchased separately from OriGene.



Anti-F11R mouse monoclonal antibody (TA506034) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY F11R ([RC200004]).