

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA506196

MSI2 Mouse Monoclonal Antibody [Clone ID: OTI2F10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2F10
Applications: IF, IHC, WB

Recommend Dilution: WB: 1:200 - 1:1000, IHC 1:150, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MSI2(NP_620412) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 35 kDa


Gene Name: musashi RNA binding protein 2

Database Link: NP 620412 Entrez Gene 124540 Human


Synonyms: MSI2H



Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MSI2 ([RC201003], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MSI2. Positive lysates [LY403371] (100ug) and [LC403371] (20ug) can be purchased separately from OriGene.




Western blot analysis of extracts (35ug) from 9 different cell lines by using anti-MSI2 monoclonal antibody (HepG2: human; HeLa: human; SVT2: mouse; A549: human; COS7: monkey; Jurkat: human; MDCK: canine; PC12: rat; MCF7: human) (1:200).






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-MSI2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506196)




Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-MSI2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506196)




Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-MSI2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506196)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-MSI2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506196)




Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-MSI2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506196)



Anti-MSI2 mouse monoclonal antibody (TA506196) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MSI2 ([RC201003]).





Immunofluorescent staining of HeLa cells using anti-MSI2 mouse monoclonal antibody (TA506196).