

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA506658**

## PMS2 Mouse Monoclonal Antibody [Clone ID: OTI2G5]

**Product data:** 

**Product Type: Primary Antibodies** 

Clone Name: OTI2G5

**Applications:** IF, IHC, WB

**Recommend Dilution:** WB 1:500, IHC 1:150, IF: 1:100

Reactivity: Human Mouse Host: Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human PMS2(NP\_000526) produced in HEK293T

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 95.6 kDa

Gene Name: PMS1 homolog 2, mismatch repair system component

NP 000526 Entrez Gene 5395 Human Database Link:

Background: This gene is one of the PMS2 gene family members found in clusters on chromosome 7. The

> product of this gene is involved in DNA mismatch repair. It forms a heterodimer with MLH1 and this complex interacts with other complexes bound to mismatched bases. Mutations in this gene are associated with hereditary nonpolyposis colorectal cancer, Turcot syndrome, and are a cause of supratentorial primitive neuroectodermal tumors. Alternatively spliced

transcript variants have been observed for this gene. [provided by RefSeq, Jul 2008]


Synonyms: HNPCC4; PMS2CL; PMSL2


**Protein Families:** Druggable Genome **Protein Pathways:** Mismatch repair






## **Product images:**









HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PMS2 ([RC210880], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PMS2. Positive lysates [LY424658] (100ug) and [LC424658] (20ug) can be purchased separately from OriGene.


Equivalent amounts of cell lysates (10 ug per lane) of wild-type HeLa cells (WT, Cat# LC810HELA) and PMS2-Knockout HeLa cells (KO, Cat# [LC831286]) were separated by SDS-PAGE and immunoblotted with anti-PMS2 monoclonal antibody TA506658 (1:500`). Then the blotted membrane was stripped and reprobed with anti-HSP90 antibody as a loading control.

Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human breast tissue using anti-PMS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506658)











Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-PMS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506658)




Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-PMS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506658)






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-PMS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506658)



Immunohistochemical staining of paraffinembedded Carcinoma of Human pancreas tissue using anti-PMS2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 120°C for 3min, TA506658)



Immunofluorescent staining of HeLa cells using anti-PMS2 mouse monoclonal antibody (TA506658).