

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA507002

MFAP3 Mouse Monoclonal Antibody [Clone ID: OTI1D3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1D3
Applications: IF, WB

Recommend Dilution: WB 1:4000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MFAP3(NP_005918) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 40 kDa

Gene Name: microfibril associated protein 3


Database Link: NP 005918 Entrez Gene 4238 Human

Synonyms: DKFZp586F0219
Protein Families: Transmembrane






Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MFAP3 ([RC204656], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MFAP3. Positive lysates [LY416962] (100ug) and [LC416962] (20ug) can be purchased separately from OriGene.



Anti-MFAP3 mouse monoclonal antibody (TA507002) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY MFAP3 ([RC204656]).