

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA507004

BORCS5 Mouse Monoclonal Antibody [Clone ID: OTI1E1]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1E1
Applications: IF, WB

Recommend Dilution: WB 1:4000, IF 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human LOH12CR1(NP_477517) produced in

HEK293T cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 22 kDa


Gene Name: BLOC-1 related complex subunit 5

Database Link: NP 477517 Entrez Gene 118426 Human


Synonyms: LOH1CR12



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY LOH12CR1 ([RC203844], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-LOH12CR1. Positive lysates [LY409223] (100ug) and [LC409223] (20ug) can be purchased separately from OriGene.



Anti-LOH12CR1 mouse monoclonal antibody (TA507004) immunofluorescent staining of COS7 cells transiently transfected by pCMV6-ENTRY LOH12CR1 ([RC203844]).