

## Product datasheet for **TA507229**


### **HORMAD2 Mouse Monoclonal Antibody [Clone ID: OTI3E4]**

#### **Product data:**

<b>Product Type:</b>	Primary Antibodies
<b>Clone Name:</b>	OTI3E4
<b>Applications:</b>	WB
<b>Recommend Dilution:</b>	WB 1:4000
<b>Reactivity:</b>	Human
<b>Host:</b>	Mouse
<b>Isotype:</b>	IgG1
<b>Clonality:</b>	Monoclonal
<b>Immunogen:</b>	Full length human recombinant protein of human HORMAD2(NP_689723) produced in HEK293T cell.
<b>Formulation:</b>	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
<b>Concentration:</b>	1 mg/ml
<b>Purification:</b>	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
<b>Predicted Protein Size:</b>	35.1 kDa
<b>Gene Name:</b>	HORMA domain containing 2
<b>Database Link:</b>	<a href="#">NP_689723 Entrez Gene 150280 Human</a>
<b>Synonyms:</b>	CT46.2



[View online »](#)

**Product images:**

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HORMAD2 ([RC205453], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HORMAD2. Positive lysates [LY407512] (100ug) and [LC407512] (20ug) can be purchased separately from OriGene.