

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA800156

TIGD1 Mouse Monoclonal Antibody [Clone ID: OTI4G9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4G9
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-267 of human

TIGD1(NP_663748) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 67.1 kDa

Gene Name: tigger transposable element derived 1

Database Link: NP 663748 Entrez Gene 200765 Human

Background: The protein encoded by this gene belongs to the tigger subfamily of the pogo superfamily of

DNA-mediated transposons in humans. These proteins are related to DNA transposons found in fungi and nematodes, and more distantly to the Tc1 and mariner transposases. They are also very similar to the major mammalian centromere protein B. The exact function of

this gene is not known. [provided by RefSeq]


Synonyms: EEYORE

Protein Families: Transcription Factors







Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TIGD1 ([RC207323], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TIGD1. Positive lysates [LY407903] (100ug) and [LC407903] (20ug) can be purchased separately from OriGene.





Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-TIGD1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA800156)




Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-TIGD1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA800156)









Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-TIGD1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 10mM citric buffer, pH6.0, 100°C for 10min, TA800156)