

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA800918**

## MAGEA3 Mouse Monoclonal Antibody [Clone ID: OTI2E11]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI2E11

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse

**Isotype:** IgG2b

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MAGEA3 (NP\_005353) produced in

HEK293T cell.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 34.6 kDa

**Gene Name:** MAGE family member A3

Database Link: NP 005353 Entrez Gene 4102 Human

**Background:** This gene is a member of the MAGEA gene family. The members of this family encode


proteins with 50 to 80% sequence identity to each other. The promoters and first exons of the MAGEA genes show considerable variability, suggesting that the existence of this gene family enables the same function to be expressed under different transcriptional controls. The MAGEA genes are clustered at chromosomal location Xq28. They have been implicated in some hereditary disorders, such as dyskeratosis congenita. [provided by RefSeq, Jul 2008]

Synonyms: CT1.3; HIP8; HYPD; MAGE3; MAGEA6





## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MAGEA3 ([RC203288], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MAGEA3. Positive lysates [LY417358] (100ug) and [LC417358] (20ug) can be purchased separately from OriGene.