

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA802185

OXSM Mouse Monoclonal Antibody [Clone ID: OTI4E10]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4E10
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG2b

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 78-343 of human OXSM

(NP_060367) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 48.7 kDa

Gene Name: 3-oxoacyl-ACP synthase, mitochondrial

Database Link: NP 060367 Entrez Gene 54995 Human

Background: This gene encodes a beta-ketoacyl synthetase. The encoded enzyme is required for

elongation of fatty acid chains in the mitochondria. Alternatively spliced transcript variants

have been described. [provided by RefSeq, Feb 2009]

Synonyms: FASN2D; KASI; KS


Protein Families: Druggable Genome

Protein Pathways: Fatty acid biosynthesis, Metabolic pathways






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY OXSM ([RC201500], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-OXSM. Positive lysates [LY413479] (100ug) and [LC413479] (20ug) can be purchased separately from OriGene.




Western blot analysis of HT29 cell lysate (35ug) by using anti-OXSM monoclonal antibody. Dilution: 1:500






Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-OXSM mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH9.0, 120°C for 3min, TA802185)



Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-OXSM mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH9.0, 120°C for 3min, TA802185)



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human endometrium tissue using anti-OXSM mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH9.0, 120°C for 3min, TA802185)