

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA802233**

## **DDIT3 Mouse Monoclonal Antibody [Clone ID: OTI1E7]**

## **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1E7
Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human DDIT3 (NP\_004074) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 19 kDa

Gene Name: DNA damage inducible transcript 3

Database Link: NP 004074 Entrez Gene 1649 Human

Background: This gene encodes a member of the CCAAT/enhancer-binding protein (C/EBP) family of

transcription factors. The protein functions as a dominant-negative inhibitor by forming heterodimers with other C/EBP members, such as C/EBP and LAP (liver activator protein), and

preventing their DNA binding activity. The protein is implicated in adipogenesis and

erythropoiesis, is activated by endoplasmic reticulum stress, and promotes apoptosis. Fusion

of this gene and FUS on chromosome 16 or EWSR1 on chromosome 22 induced by translocation generates chimeric proteins in myxoid liposarcomas or Ewing sarcoma.

Multiple alternatively spliced transcript variants encoding two isoforms with different length

have been identified. [provided by RefSeq, Aug 2010]

Synonyms: CEBPZ; CHOP; CHOP-10; CHOP10; GADD153


Protein Families: Druggable Genome, Transcription Factors

**Protein Pathways:** MAPK signaling pathway





# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DDIT3 ([RC201301], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DDIT3. Positive lysates [LY418226] (100ug) and [LC418226] (20ug) can be purchased separately from OriGene.