

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA802368

SPATA2L Mouse Monoclonal Antibody [Clone ID: OTI3F6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI3F6

Applications: IHC, WB

Recommend Dilution: WB 1:2000

Reactivity: Human

Host: Mouse IgG2a

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 179-424 of human

SPATA2L (NP_689552) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 46 kDa

Gene Name: spermatogenesis associated 2 like


Database Link: NP 689552 Entrez Gene 124044 Human

Synonyms: C16orf76; tamo






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SPATA2L ([RC208075], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SPATA2L. Positive lysates [LY407590] (100ug) and [LC407590] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-SPATA2L mouse monoclonal antibody. (TA802368) Dilution: 1:150






Immunohistochemical staining of paraffinembedded Carcinoma of Human thyroid tissue using anti-SPATA2L mouse monoclonal antibody. (TA802368) Dilution: 1:150



Immunohistochemical staining of paraffinembedded Human lymph node tissue within the normal limits using anti-SPATA2L mouse monoclonal antibody. (TA802368) Dilution: 1:150



Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-SPATA2L mouse monoclonal antibody. (TA802368) Dilution: 1:150