

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA802721**

## PARN Mouse Monoclonal Antibody [Clone ID: OTI2D9]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI2D9

Applications: IHC, WB

Recommend Dilution: WB 1:2000

Reactivity: Human

Host: Mouse

Isotype: IgG2b

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 1-266 of human PARN

(NP\_002573) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 73.3 kDa

**Gene Name:** poly(A)-specific ribonuclease

Database Link: NP 002573 Entrez Gene 5073 Human

**Background:** The protein encoded by this gene is a 3'-exoribonuclease, with similarity to the RNase D

family of 3'-exonucleases. It prefers poly(A) as the substrate, hence, efficiently degrades poly(A) tails of mRNAs. Exonucleolytic degradation of the poly(A) tail is often the first step in the decay of eukaryotic mRNAs. This protein is also involved in silencing of certain maternal mRNAs during oocyte maturation and early embryonic development, as well as in nonsensemediated decay (NMD) of mRNAs that contain premature stop codons. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by

RefSeq, Aug 2008]

Synonyms: DAN


**Protein Families:** Transcription Factors

**Protein Pathways:** RNA degradation







## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY PARN ([RC207220], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PARN. Positive lysates [LY419233] (100ug) and [LC419233] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human ovary tissue using anti-PARN mouse monoclonal antibody. (TA802721) Dilution: 1:150



Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-PARN mouse monoclonal antibody. (TA802721) Dilution: 1:150