

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA803200

ROR alpha (RORA) Mouse Monoclonal Antibody [Clone ID: OTI2E8]

Product data:

Clone Name:OTIZE8Applications:WBRecommed Dilution:WB 1:500Reactivity:HumanHost:MouseIsotype:IgG2aClonality:MonoclonalImmunogen:Human recombinant protein fragment corresponding to amino acids 98-345 of human RORA (NP_002934) produced in E.coli.Formulation:PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.Oncentration:Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein SizeAR related orphan receptor ADatabase Link:Ne pro2334 Entrez Gene 6095 Human variant of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside shiphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. Four transcript variants encoding different isoforms have been described for this gene. Four transcript variants encoding different isoforms have been described for this gene. Four transcript variants encoding different isoforms have been described for this gene. Four transcript variants encoding different isoforms have been described for this gene. Four transcript variants encoding different isoforms have been described for this gene. Four transcript variants encoding different isoforms have been described for this gene. Four transcript variants encoding different isoforms have been described for this gene. Four transcript variants encoding different isoforms have been described for this gene.	Product Type:	Primary Antibodies
Recommend Dilution:WB 1:500Reactivity:HumanHost:MouseIsotype:IgG2aClonality:MonoclonalImmunogen:Human recombinant protein fragment corresponding to amino acids 98-345 of human RORA (NP_002934) produced in E.coli.Formulation:PS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.Concentration:Img/mlPurified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein Size:62.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq. Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Clone Name:	OTI2E8
Reactivity:HumanHost:MouseIsotype:IgG2aClonality:MonoclonalImmunogen:Human recombinant protein fragment corresponding to amino acids 98-345 of human RORA (NP_002934) produced in E.coli.Formulation:PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.Forneutation:1 mg/mlPurification:Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein Size6.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those genes. The specific functions upstream of several genes to enhance the expression of those g	Applications:	WB
Host:MouseIsotype:IgG2aClonality:MonoclonalImmunogen:Human recombinant protein fragment corresponding to amino acids 98-345 of human RORA (NP_002934) produced in E.coli.Formulation:PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.Concentration:1 mg/mlPurification:Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein Size:62.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Recommend Dilution:	WB 1:500
Isotype:IgG2aClonality:MonoclonalImmunogen:Human recombinant protein fragment corresponding to amino acids 98-345 of human RORA (NP_002934) produced in E.coli.Formulation:PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.Concentration:1 mg/mlPurification:Quified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein Size:62.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:Ne protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq. Jul 2008]Synonyms:NR1E1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Reactivity:	Human
Clonality:MonoclonalImmunogen:Human recombinant protein fragment corresponding to amino acids 98-345 of human RORA (NP_002934) produced in E.coli.Formulation:PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.Concentration:1 mg/mlPurification:Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein Size:62.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Host:	Mouse
Immunogen:Human recombinant protein fragment corresponding to amino acids 98-345 of human RORA (NP_002934) produced in E.coli.Formulation:PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.Concentration:1 mg/mlPurification:Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein Size:62.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	lsotype:	lgG2a
(NP_002934) produced in E.coli.Formulation:PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.Concentration:1 mg/mlPurification:Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein Size:62.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP_002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Clonality:	Monoclonal
Concentration:1 mg/mlPurification:Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein Size:62.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Immunogen:	
Purification:Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)Predicted Protein Size:62.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
(protein A/G)Predicted Protein Size:62.1 kDaGene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Concentration:	1 mg/ml
Gene Name:RAR related orphan receptor ADatabase Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Purification:	
Database Link:NP 002934 Entrez Gene 6095 HumanBackground:The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008]Synonyms:NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA	Predicted Protein Size:	62.1 kDa
 Background: The protein encoded by this gene is a member of the NR1 subfamily of nuclear hormone receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008] Synonyms: NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA 	Gene Name:	RAR related orphan receptor A
 receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008] Synonyms: NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA 	Database Link:	<u>NP_002934 Entrez Gene 6095 Human</u>
	Background:	receptors. It can bind as a monomer or as a homodimer to hormone response elements upstream of several genes to enhance the expression of those genes. The specific functions of this protein are not known, but it has been shown to interact with NM23-2, a nucleoside diphosphate kinase involved in organogenesis and differentiation, as well as with NM23-1, the product of a tumor metastasis suppressor candidate gene. Four transcript variants
Protein Families: Druggable Genome, Nuclear Hormone Receptor, Transcription Factors	Synonyms:	NR1F1; ROR1; ROR2; ROR3; RZR-ALPHA; RZRA
	Protein Families:	Druggable Genome, Nuclear Hormone Receptor, Transcription Factors



View online »

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



Product images:

_	
_	
_	
	-
—	
—	
—	
_	

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY RORA ([RC214618], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-RORA. Positive lysates [LY401041] (100ug) and [LC401041] (20ug) can be purchased separately from OriGene.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US