

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# Product datasheet for TA803430

## Suppressor of Ty 4 homolog 1 (SUPT4H1) Mouse Monoclonal Antibody [Clone ID: OTI6B10]

#### **Product data:**

Product Type:	Primary Antibodies
Clone Name:	OTI6B10
Applications:	WB
<b>Recommend Dilution:</b>	WB 1:500
Reactivity:	Human
Host:	Mouse
lsotype:	lgG1
Clonality:	Monoclonal
Immunogen:	Full length human recombinant protein of human SUPT4H1 (NP_003159) produced in E.coli.
Formulation:	PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.
Concentration:	1 mg/ml
Purification:	Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography (protein A/G)
Predicted Protein Size:	13 kDa
Gene Name:	SPT4 homolog, DSIF elongation factor subunit
Database Link:	<u>NP_003159 Entrez Gene 6827 Human</u>
Synonyms:	SPT4; SPT4H; Supt4a; SUPT4H
Protein Families:	Transcription Factors



This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US



## **Product images:**

170	-	
130	-	
100	-	
70	-	
55	-	
40		
35	-	
25	-	
15	-	-
10	-	

HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SUPT4H1 ([RC201211], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SUPT4H1. Positive lysates [LY418857] (100ug) and [LC418857] (20ug) can be purchased separately from OriGene.

This product is to be used for laboratory only. Not for diagnostic or therapeutic use. ©2020 OriGene Technologies, Inc., 9620 Medical Center Drive, Ste 200, Rockville, MD 20850, US