

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA803473

GBA Mouse Monoclonal Antibody [Clone ID: OTI2C12]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2C12
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:100

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 40-315 of human GBA

(NP_000148) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 55.5 kDa

Gene Name: glucosylceramidase beta

Database Link: NP 000148 Entrez Gene 2629 Human

Background: This gene encodes a lysosomal membrane protein that cleaves the beta-glucosidic linkage of

glycosylceramide, an intermediate in glycolipid metabolism. Mutations in this gene cause

Gaucher disease, a lysosomal storage disease characterized by an accumulation of

glucocerebrosides. A related pseudogene is approximately 12 kb downstream of this gene on

chromosome 1. Alternative splicing results in multiple transcript variants. [provided by

RefSeq, Jan 2010]

Synonyms: GBA1; GCB; GLUC


Protein Families: Druggable Genome

Protein Pathways: Lysosome, Metabolic pathways, Other glycan degradation, Sphingolipid metabolism






Product images:




E.coli lysate (left lane) and E.coli lysate expressing human recombinant protein fragment corresponding to amino acids 40-315 of human GBA (NP_000148) were separated by SDS-PAGE and immunoblotted with anti-GBA.




Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-GBA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA803473)






Immunohistochemical staining of paraffinembedded Human pancreas tissue within the normal limits using anti-GBA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA803473)



Immunohistochemical staining of paraffinembedded Human thyroid tissue within the normal limits using anti-GBA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA803473)



Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-GBA mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA803473)