

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA803844**

### Angiotensin Converting Enzyme 2 (ACE2) Mouse Monoclonal Antibody [Clone ID: OTI4C5]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4C5
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 18-237 of human ACE2

(NP\_068576) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Gene Name:** angiotensin I converting enzyme 2

Database Link: NP 068576 Entrez Gene 59272 Human

**Background:** The protein encoded by this gene belongs to the angiotensin-converting enzyme family of

dipeptidyl carboxydipeptidases and has considerable homology to human angiotensin 1 converting enzyme. This secreted protein catalyzes the cleavage of angiotensin I into angiotensin 1-9, and angiotensin II into the vasodilator angiotensin 1-7. The organ- and cell-

specific expression of this gene suggests that it may play a role in the regulation of

cardiovascular and renal function, as well as fertility. In addition, the encoded protein is a functional receptor for the spike glycoprotein of the human coronaviruses SARS and HCoV-

NL63. [provided by RefSeq, Jul 2008]

Synonyms: ACEH


**Protein Families:** Druggable Genome, Secreted Protein, Transmembrane

**Protein Pathways:** Renin-angiotensin system






## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ACE2 ([RC208442], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ACE2.



Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-ACE2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA803844)