

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA804063

L Kynurenine Hydrolase (KYNU) Mouse Monoclonal Antibody [Clone ID: OTI5E4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI5E4
Applications: IHC, WB

Recommend Dilution: WB 1:500, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-216 of human KYNU

(NP_001028170) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 34.5 kDa

Gene Name: kynureninase

Database Link: NP 001028170 Entrez Gene 8942 Human

Background: Kynureninase is a pyridoxal-5'-phosphate (pyridoxal-P) dependent enzyme that catalyzes the

cleavage of L-kynurenine and L-3-hydroxykynurenine into anthranilic and 3-

hydroxyanthranilic acids, respectively. Kynureninase is involved in the biosynthesis of NAD cofactors from tryptophan through the kynurenine pathway. Alternative splicing results in

multiple transcript variants. [provided by RefSeq, Nov 2010]

Synonyms: KYNUU


Protein Families: Protease

Protein Pathways: Metabolic pathways, Tryptophan metabolism






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY KYNU ([RC201559], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-KYNU. Positive lysates [LY422334] (100ug) and [LC422334] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Human breast tissue within the normal limits using anti-KYNU mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA804063)





Immunohistochemical staining of paraffinembedded Carcinoma of Human liver tissue using anti-KYNU mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA804063)



Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-KYNU mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA804063)