

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA804627

EPOR Mouse Monoclonal Antibody [Clone ID: OTI4F12]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4F12
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 25-250 of human EPOR

(NP_000112) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Gene Name: erythropoietin receptor

Database Link: NP 000112 Entrez Gene 2057 Human

Background: This gene encodes the erythropoietin receptor which is a member of the cytokine receptor

family. Upon erythropoietin binding, this receptor activates Jak2 tyrosine kinase which activates different intracellular pathways including: Ras/MAP kinase, phosphatidylinositol 3-kinase and STAT transcription factors. The stimulated erythropoietin receptor appears to have a role in erythroid cell survival. Defects in the erythropoietin receptor may produce erythroleukemia and familial erythrocytosis. Dysregulation of this gene may affect the growth of certain tumors. Alternate splicing results in multiple transcript variants. [provided by

RefSeq, May 2010]

Synonyms: EPO-R

Protein Families: Druggable Genome, Secreted Protein, Transmembrane


Protein Pathways: Cytokine-cytokine receptor interaction, Hematopoietic cell lineage, Jak-STAT signaling

pathway






Product images:



Human recombinant protein fragment corresponding to amino acids 25-250 of human EPOR (NP_000112) produced in E.coli.



Immunohistochemical staining of paraffinembedded Human endometrium tissue within the normal limits using anti-EPOR mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA804627)