

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA805346

HIP55 (DBNL) Mouse Monoclonal Antibody [Clone ID: OTI7B7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI7B7
Applications: WB

December of Biletians N/D 1:200

Recommend Dilution: WB 1:2000
Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human DBNL (NP_054782) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)


Predicted Protein Size: 48.1 kDa **Gene Name:** drebrin like

Database Link: NP 054782 Entrez Gene 28988 Human

Synonyms: ABP1; HIP-55; HIP55; SH3P7



Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY DBNL ([RC203435], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-DBNL. Positive lysates [LY402274] (100ug) and [LC402274] (20ug) can be purchased separately from OriGene.