

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA805784

MCM2 Mouse Monoclonal Antibody [Clone ID: OTI1B1]

Product data:

Isotype:

Product Type: Primary Antibodies

Clone Name: OTI1B1
Applications: IHC

Recommend Dilution: IHC 1:150
Reactivity: Human
Host: Mouse

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human MCM2 (NP_004517) produced in HEK293T

cell.

lgG1

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 101.7 kDa

Gene Name: minichromosome maintenance complex component 2

Database Link: NP 004517 Entrez Gene 4171 Human

Background: The protein encoded by this gene is one of the highly conserved mini-chromosome

maintenance proteins (MCM) that are involved in the initiation of eukaryotic genome replication. The hexameric protein complex formed by MCM proteins is a key component of the pre-replication complex (pre_RC) and may be involved in the formation of replication forks and in the recruitment of other DNA replication related proteins. This protein forms a complex with MCM4, 6, and 7, and has been shown to regulate the helicase activity of the complex. This protein is phosphorylated, and thus regulated by, protein kinases CDC2 and CDC7. Multiple alternatively spliced transcript variants have been found, but the full-length

nature of some variants has not been defined. [provided by RefSeq, Oct 2012]

Synonyms: BM28; CCNL1; cdc19; CDCL1; D3S3194; MITOTIN


Protein Families: Druggable Genome, Stem cell - Pluripotency, Transcription Factors






Protein Pathways: Cell cycle, DNA replication

Product images:




Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-MCM2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA805784)



Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-MCM2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA805784)





Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-MCM2 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA805784)