

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA805797

Interferon gamma (IFNG) Mouse Monoclonal Antibody [Clone ID: OTI1G7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI1G7

Applications: WB

Recommend Dilution: WB 1:2000 Reactivity: Human

Mouse Host:

Isotype: lgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 24-166 of human

IFNG(NP 000610) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 16.7 kDa

Gene Name: interferon, gamma

Database Link: NP 000610 Entrez Gene 3458 Human

Background: This gene encodes a member of the type II interferon family. The protein encoded is a

soluble cytokine with antiviral, immunoregulatory and anti-tumor properties and is a potent

activator of macrophages. Mutations in this gene are associated with aplastic anemia.

[provided by RefSeq, Nov 2009]

Synonyms: IFG: IFI

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: Allograft rejection, Cytokine-cytokine receptor interaction, Graft-versus-host disease, Jak-STAT


> signaling pathway, Natural killer cell mediated cytotoxicity, Proteasome, Regulation of autophagy, Systemic lupus erythematosus, T cell receptor signaling pathway, TGF-beta

signaling pathway, Type I diabetes mellitus





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY IFNG ([RC209993], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-IFNG. Positive lysates [LY400207] (100ug) and [LC400207] (20ug) can be purchased separately from OriGene.