

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA805877

hSET1 (SETD1A) Mouse Monoclonal Antibody [Clone ID: OTI7B7]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI7B7
Applications: IHC, WB

Recommend Dilution: WB 1:500, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 1-320 of human

SETD1A(NP_055527) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 185.9 kDa

Gene Name: SET domain containing 1A

Database Link: NP 055527 Entrez Gene 9739 Human

Background: SET1A is a component of a histone methyltransferase (HMT) complex that produces mono-,

di-, and trimethylated histone H3 at Lys4. The complex is the analog of the S. cerevisiae Set1/COMPASS complex (Lee and Skalnik, 2005 [PubMed 16253997]). Also see SET1B (MIM

611055). [supplied by OMIM, Mar 2008]

Synonyms: KMT2F; Set1; Set1A


Protein Families: Druggable Genome

Protein Pathways: Lysine degradation






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY SETD1A ([RC214996], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-SETD1A. Positive lysates [LY415047] (100ug) and [LC415047] (20ug) can be purchased separately from OriGene.



Immunohistochemical staining of paraffinembedded Human prostate tissue within the normal limits using anti-SETD1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA805877)





Immunohistochemical staining of paraffinembedded Human tonsil within the normal limits using anti-SETD1A mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA805877)