

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA806417**

## Lipin 1 (LPIN1) Mouse Monoclonal Antibody [Clone ID: OTI4C11]

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4C11

Applications: IHC

Recommend Dilution: IHC 1:150

Reactivity: Human Host: Mouse

**Isotype:** lgG2b

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 1-274 of human

LPIN1(NP\_663731) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 98.5 kDa **Gene Name:** lipin 1

Database Link: NP 663731 Entrez Gene 23175 Human

**Background:** This gene encodes a magnesium-ion-dependent phosphatidic acid phosphohydrolase

enzyme that catalyzes the penultimate step in triglyceride synthesis including the dephosphorylation of phosphatidic acid to yield diacylglycerol. Expression of this gene is required for adipocyte differentiation and it also functions as a nuclear transcriptional

coactivator with some peroxisome proliferator-activated receptors to modulate expression of

other genes involved in lipid metabolism. Mutations in this gene are associated with metabolic syndrome, type 2 diabetes, and autosomal recessive acute recurrent

myoglobinuria (ARARM). This gene is also a candidate for several human lipodystrophy syndromes. Alternative splicing results in multiple transcript variants encoding distinct isoforms. Additional splice variants have been described but their full-length structures have


not been determined. [provided by RefSeq, May 2012]






Synonyms: PAP1

# **Product images:**




Immunohistochemical staining of paraffinembedded Human breast tissue within the normal limits using anti-LPIN1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA806417)




Immunohistochemical staining of paraffinembedded Human colon tissue within the normal limits using anti-LPIN1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA806417)





Immunohistochemical staining of paraffinembedded Human Kidney tissue within the normal limits using anti-LPIN1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA806417)



Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-LPIN1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1 mM EDTA in 10mM Tris, pH8.5, 120°C for 3min, TA806417)