

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA806458

TNFRSF18 Mouse Monoclonal Antibody [Clone ID: OTI4A9]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4A9
Applications: FC, WB

Recommend Dilution: WB 1:2000

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 26-162 of human

TNFRSF18(NP_004186) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 23.5 kDa

Gene Name: tumor necrosis factor receptor superfamily member 18

Database Link: NP 004186 Entrez Gene 8784 Human

Background: This gene encodes a member of the TNF-receptor superfamily. The encoded receptor has

been shown to have increased expression upon T-cell activation, and it is thought to play a key role in dominant immunological self-tolerance maintained by CD25(+)CD4(+) regulatory T cells. Knockout studies in mice also suggest the role of this receptor is in the regulation of CD3-driven T-cell activation and programmed cell death. Three alternatively spliced transcript variants of this gene encoding distinct isoforms have been reported. [provided by RefSeq, Feb

Synonyms: AITR; CD357; GITR; GITR-D


Protein Families: Druggable Genome, Secreted Protein, Transmembrane

Protein Pathways: Cytokine-cytokine receptor interaction







Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY TNFRSF18 ([RC219997], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-TNFRSF18. Positive lysates [LY418153] (100ug) and [LC418153] (20ug) can be purchased separately from OriGene.



Western blot analysis of human peripheral blood mononuclear cells (PBMCs) by using anti-TNFRSF18 monoclonal antibody. Left: Unstimulated PBMC. Right: PBMC stimulated with 10ug/ml Phytohemagglutinin for 3 days.





Flow cytometric analysis of TNFRSF18 expression on living stimulated human peripheral blood mononuclear cells (PBMCs). Unstimulated PBMCs, (a.left, b.blue), or 10ug/ml Phytohemagglutinin-stimulated (3 days) PBMCs, (a.right, b.red), were stained with Anti-Human TNFRSF18 antibody (TA806458) (1:100).