

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA806787**

## NR2F6 Mouse Monoclonal Antibody [Clone ID: OTI2H4]

### **Product data:**

Isotype:

**Product Type:** Primary Antibodies

Clone Name: OTI2H4

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human

**Host:** Mouse

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human NR2F6 (NP\_005225) produced in HEK293T

cell.

IgG2a

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 42.8 kDa

**Gene Name:** nuclear receptor subfamily 2 group F member 6


Database Link: NP 005225 Entrez Gene 2063 Human

**Synonyms:** EAR-2; EAR2; ERBAL2

**Protein Families:** Druggable Genome, Nuclear Hormone Receptor, Transcription Factors



# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY NR2F6 ([RC201336], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-NR2F6. Positive lysates [LY401603] (100ug) and [LC401603] (20ug) can be purchased separately from OriGene.