

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA806909

Integrin beta 1 (ITGB1) Mouse Monoclonal Antibody [Clone ID: OTI2A3]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2A3

Applications: FC, WB

Recommend Dilution: WB 1:500

Reactivity: Human

Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ITGB1 (NP_002202) produced in HEK293T

cell

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 86.1 kDa

Gene Name: integrin subunit beta 1

Database Link: NP 002202 Entrez Gene 3688 Human

Synonyms: CD29; FNRB; GPIIA; MDF2; MSK12; VLA-BETA; VLAB

Protein Families: Druggable Genome, ES Cell Differentiation/IPS, Transmembrane

Protein Pathways: Arrhythmogenic right ventricular cardiomyopathy (ARVC), Axon guidance, Cell adhesion

molecules (CAMs), Dilated cardiomyopathy, ECM-receptor interaction, Focal adhesion, Hypertrophic cardiomyopathy (HCM), Leukocyte transendothelial migration, Pathogenic


Escherichia coli infection, Pathways in cancer, Regulation of actin cytoskeleton, Small cell lung

cancer






Product images:





HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ITGB1 ([RC203818], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ITGB1 (1:500).




Western blot analysis of extracts (35ug) from 3 different cell lines by using anti-ITGB1 monoclonal antibody (1:500).











Flow cytometric Analysis of living Hela cells, using anti-ITGB1 antibody (TA806909), (Red), compared to an IgG isotype control, (green), and negative control (PBS), (Blue) (1:100).



Flow cytometric Analysis of living U-87MG cells, using anti-ITGB1 antibody (TA806909), (Red), compared to an IgG isotype control, (green), and negative control (PBS), (Blue) (1:100).





Flow cytometric Analysis of living 786-O cells, using anti-ITGB1 antibody (TA806909), (Red), compared to an IgG isotype control, (green), and negative control (PBS), (Blue) (1:100).