

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA807045

Adrenomedullin (ADM) Mouse Monoclonal Antibody [Clone ID: OTI2A4]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2A4

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 22-146 of human

ADM(NP_001115) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 20.2 kDa

Gene Name: adrenomedullin

Database Link: NP 001115 Entrez Gene 133 Human

Background: Adrenomedullin, a hypotensive peptide found in human pheochromocytoma, consists of 52

amino acids, has 1 intramolecular disulfide bond, and shows a slight homology with the calcitonin gene-related peptide. It may function as a hormone in circulation control because it

is found in blood in a considerable concentration. The precursor, called

preproadrenomedullin, is 185 amino acids long. By RNA-blot analysis, human adrenomedullin mRNA was found to be highly expressed in several tissues. Genomic ADM DNA consists of 4 exons and 3 introns, with the 5-prime flanking region containing TATA, CAAT, and GC boxes. There are also multiple binding sites for activator protein-2 and a cAMP-regulated enhancer

element. [provided by RefSeq, Jul 2008]


Synonyms: AM; PAMP

Protein Families: Druggable Genome, Secreted Protein





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ADM ([RC200507], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ADM. Positive lysates [LY420111] (100ug) and [LC420111] (20ug) can be purchased separately from OriGene.