

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA807260

HIPK1 Mouse Monoclonal Antibody [Clone ID: OTI2C11]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2C11

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

Immunogen: Human recombinant protein fragment corresponding to amino acids 439-697 of human

HIPK1(NP_852003) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 87.1 kDa

Gene Name: homeodomain interacting protein kinase 1

Database Link: NP 852003 Entrez Gene 204851 Human

Background: The protein encoded by this gene belongs to the Ser/Thr family of protein kinases and HIPK

subfamily. It phosphorylates homeodomain transcription factors and may also function as a co-repressor for homeodomain transcription factors. Alternative splicing results in four transcript variants encoding four distinct isoforms. [provided by RefSeq, Jul 2008]


Synonyms: Myak; Nbak2

Protein Families: Druggable Genome, Protein Kinase, Transcription Factors





Product images:



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY HIPK1 ([RC206176], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-HIPK1 (1:2000). Positive lysates [LY405750] (100ug) and [LC405750] (20ug) can be purchased separately from OriGene.