

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA807444**

## LOXL2 Mouse Monoclonal Antibody [Clone ID: OTI4H3]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI4H3

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

**Isotype:** IgG2b

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 152-450 of human

LOXL2(NP\_002309) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 84 kDa

**Gene Name:** lysyl oxidase like 2

Database Link: NP 002309 Entrez Gene 4017 Human

**Background:** This gene encodes a member of the lysyl oxidase gene family. The prototypic member of the

family is essential to the biogenesis of connective tissue, encoding an extracellular copperdependent amine oxidase that catalyses the first step in the formation of crosslinks in collagens and elastin. A highly conserved amino acid sequence at the C-terminus end appears to be sufficient for amine oxidase activity, suggesting that each family member may retain this function. The N-terminus is poorly conserved and may impart additional roles in

developmental regulation, senescence, tumor suppression, cell growth control, and

chemotaxis to each member of the family. [provided by RefSeq, Jul 2008]


Synonyms: LOR2; WS9-14


**Protein Families:** Druggable Genome, Secreted Protein





### **Product images:**





HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY LOXL2 ([RC200455], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-LOXL2 (1:2000). Positive lysates [LY400844] (100ug) and [LC400844] (20ug) can be purchased separately from OriGene.

Equivalent amounts of cell lysates (10 ug per lane) of wild-type Hela cells (WT, Cat# LC810HELA) and LOXL2-Knockout Hela cells (KO, Cat# [LC810132]) were separated by SDS-PAGE and immunoblotted with anti-LOXL2 monoclonal antibody TA807444. Then the blotted membrane was stripped and reprobed with anti-b-actin antibody ([TA811000]) as a loading control (1:500).