

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA807639**

### POLA2 Mouse Monoclonal Antibody [Clone ID: OTI2B4]

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI2B4

**Applications:** WB

**Recommend Dilution:** WB 1:2000

Reactivity: Human Host: Mouse

**Isotype:** IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 1-284 of human

POLA2(NP\_002680) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 65.8 kDa

**Gene Name:** polymerase (DNA) alpha 2, accessory subunit


Database Link: NP 002680 Entrez Gene 23649 Human

**Synonyms:** FLJ21662; FLJ37250

**Protein Pathways:** DNA replication, Metabolic pathways, Purine metabolism, Pyrimidine metabolism



## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY POLA2 ([RC201145], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-POLA2. Positive lysates [LY419168] (100ug) and [LC419168] (20ug) can be purchased separately from OriGene.