

## OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA807775**

## Midkine (MDK) Mouse Monoclonal Antibody [Clone ID: OTI8C6]

**Product data:** 

**Product Type:** Primary Antibodies

Clone Name: OTI8C6

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 21-143 of human

MDK(NP\_001012334) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 15.4 kDa

**Gene Name:** midkine (neurite growth-promoting factor 2)

Database Link: NP 001012334 Entrez Gene 4192 Human

**Background:** This gene encodes a member of a small family of secreted growth factors that binds heparin

and responds to retinoic acid. The encoded protein promotes cell growth, migration, and angiogenesis, in particular during tumorigenesis. This gene has been targeted as a therapeutic for a variety of different disorders. Alternatively spliced transcript variants

encoding multiple isoforms have been observed. [provided by RefSeq, Jul


**Synonyms:** ARAP; MK; NEGF2

**Protein Families:** Druggable Genome, Secreted Protein, Transmembrane





## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MDK ([RC221818], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MDK (1:2000). Positive lysates [LY425316] (100ug) and [LC425316] (20ug) can be purchased separately from OriGene.