

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA808044**

## **ZIC3 Mouse Monoclonal Antibody [Clone ID: OTI1G7]**

## **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1G7

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human

Host: Mouse Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 242-467 of human

ZIC3(NP\_003404) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 50.4 kDa

**Gene Name:** Zic family member 3

Database Link: NP 003404 Entrez Gene 7547 Human

**Background:** This gene encodes a member of the ZIC family of C2H2-type zinc finger proteins. This nuclear

protein probably functions as a transcription factor in early stages of left-right body axis formation. Mutations in this gene cause X-linked visceral heterotaxy, which includes

congenital heart disease and left-right axis defects in organs. [provided by RefSeq, Jul 2008]


Synonyms: HTX; HTX1; VACTERLX; ZNF203

**Protein Families:** Druggable Genome, Embryonic stem cells, ES Cell Differentiation/IPS





# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZIC3 ([RC220375], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZIC3 (1:2000). Positive lysates [LY418705] (100ug) and [LC418705] (20ug) can be purchased separately from OriGene.