

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA808408**

## FAM104B Mouse Monoclonal Antibody [Clone ID: OTI1D10]

#### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI1D10

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

**Isotype:** IgG1

Clonality: Monoclonal

**Immunogen:** Full length human recombinant protein of human FAM104B(NP\_612371) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Predicted Protein Size:** 12.9 kDa


**Gene Name:** family with sequence similarity 104 member B

Database Link: NP 612371 Entrez Gene 90736 Human

Synonyms: CXorf44



# **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY FAM104B ([RC205045], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-FAM104B (1:2000). Positive lysates [LY408633] (100ug) and [LC408633] (20ug) can be purchased separately from OriGene.