

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA808530

ZNF449 Mouse Monoclonal Antibody [Clone ID: OTI4C6]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI4C6
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG1

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ZNF449 (NP_689908) produced in E.coli.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 59.8 kDa

Gene Name: zinc finger protein 449

Database Link: NP 689908 Entrez Gene 203523 Human

Background: This gene encodes a nuclear protein that likely functions as a transcription factor. The

protein includes an N-terminal SCAN domain, and seven C2H2-type zinc finger motifs.

[provided by RefSeq, May 2010]


Synonyms: ZSCAN19

Protein Families: Transcription Factors






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ZNF449 ([RC209172], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ZNF449 (1:2000). Positive lysates [LY403487] (100ug) and [LC403487] (20ug) can be purchased separately from OriGene.




Immunohistochemical staining of paraffinembedded Carcinoma of Human kidney tissue using anti-ZNF449 mouse monoclonal antibody. (Heat-induced epitope retrieval by Tris-EDTA, pH8.0, TA808530) (1:150)





Immunohistochemical staining of paraffinembedded Carcinoma of Human lung tissue using anti-ZNF449 mouse monoclonal antibody. (Heat-induced epitope retrieval by Tris-EDTA, pH8.0, TA808530) (1:150)



Immunohistochemical staining of paraffinembedded Human bladder tissue within the normal limits using anti-ZNF449 mouse monoclonal antibody. (Heat-induced epitope retrieval by Tris-EDTA, pH8.0, TA808530) (1:150)