

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

## **Product datasheet for TA808743**

## KAL1 (ANOS1) Mouse Monoclonal Antibody [Clone ID: OTI5E10]

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI5E10

Applications: IHC

Reactivity: IHC 1:500

Human

Host: Mouse

**Isotype:** IgG2a

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 129-373 of human

KAL1(NP\_000207) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

**Gene Name:** anosmin 1

Database Link: NP 000207 Entrez Gene 3730 Human

**Background:** Mutations in this gene cause the X-linked Kallmann syndrome. The encoded protein is similar

in sequence to proteins known to function in neural cell adhesion and axonal migration. In addition, this cell surface protein is N-glycosylated and may have anti-protease activity.


[provided by RefSeq, Jul 2008]

Synonyms: ADMLX; HH1; HHA; KAL; KAL1; KALIG-1; KMS; WFDC19





# **Product images:**



Immunohistochemical staining of paraffinembedded Human embryonic cerebellum within the normal limits using anti-KAL1 mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA808743) (1:500)