

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Product datasheet for TA808779

ICOS Ligand (ICOSLG) Mouse Monoclonal Antibody [Clone ID: OTI2C5]

Product data:

Product Type: Primary Antibodies

Clone Name: OTI2C5
Applications: IHC, WB

Recommend Dilution: WB 1:2000, IHC 1:150

Reactivity: Human
Host: Mouse
Isotype: IgG2a

Clonality: Monoclonal

Immunogen: Full length human recombinant protein of human ICOSLG (NP_056074) produced in HEK293T

cell.

Formulation: PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

Purification: Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 33.2 kDa

Gene Name: inducible T-cell costimulator ligand

Database Link: NP 056074 Entrez Gene 23308 Human

Background: Ligand for the T-cell-specific cell surface receptor ICOS. Acts as a costimulatory signal for T-

cell proliferation and cytokine secretion; induces also B-cell proliferation and differentiation into plasma cells. Could play an important role in mediating local tissue responses to inflammatory conditions, as well as in modulating the secondary immune response by co-

stimulating memory T-cell function. [UniProtKB/Swiss-Prot Function]

Synonyms: B7-H2; B7h; B7H2; B7RP-1; B7RP1; CD275; GL50; ICOS-L; ICOSL; LICOS


Protein Families: Transmembrane

Protein Pathways: Cell adhesion molecules (CAMs)






Product images:




HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY ICOSLG ([RC208975], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-ICOSLG (1:2000). Positive lysates [LY402420] (100ug) and [LC402420] (20ug) can be purchased separately from OriGene.




Western blot analysis of 35ug lysates of Ramos, Raji, K562, Jurkat cell lines and TNF-a recombinant protein (100ng/ml, [TP750007]) sitimulated HUVEC cells for 0h, 6h and 16h (1:500).






Immunohistochemical staining of paraffinembedded Adenocarcinoma of Human colon tissue using anti-ICOSLG mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA808779) (1:150)



Immunohistochemical staining of paraffinembedded Human liver tissue within the normal limits using anti-ICOSLG mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA808779) (1:150)



Immunohistochemical staining of paraffinembedded Carcinoma of Human prostate tissue using anti-ICOSLG mouse monoclonal antibody. (Heat-induced epitope retrieval by 1mM EDTA in 10mM Tris buffer (pH8.5) at 120°C for 3min, TA808779) (1:150)