

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# **Product datasheet for TA808790**

## MAST205 (MAST2) Mouse Monoclonal Antibody [Clone ID: OTI3A5]

### **Product data:**

**Product Type:** Primary Antibodies

Clone Name: OTI3A5

Applications: WB

Recommend Dilution: WB 1:2000

Reactivity: Human Host: Mouse

Isotype: IgG1

Clonality: Monoclonal

**Immunogen:** Human recombinant protein fragment corresponding to amino acids 1482-1798 of human

MAST2(NP\_055927) produced in E.coli.

**Formulation:** PBS (PH 7.3) containing 1% BSA, 50% glycerol and 0.02% sodium azide.

Concentration: 1 mg/ml

**Purification:** Purified from mouse ascites fluids or tissue culture supernatant by affinity chromatography

(protein A/G)

Predicted Protein Size: 196.3 kDa

**Gene Name:** microtubule associated serine/threonine kinase 2


Database Link: NP 055927 Entrez Gene 23139 Human

Synonyms: MAST205; MTSSK


**Protein Families:** Druggable Genome, Protein Kinase



## **Product images:**



HEK293T cells were transfected with the pCMV6-ENTRY control (Left lane) or pCMV6-ENTRY MAST2 ([RC209492], Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-MAST2 (1:2000). Positive lysates [LY402408] (100ug) and [LC402408] (20ug) can be purchased separately from OriGene.



Equivalent amounts of cell lysates (10 ug per lane) of wild-type HeLa cells (WT, Cat# LC810HELA) and MAST2-Knockout HeLa cells (KO, Cat# [LC830654]) were separated by SDS-PAGE and immunoblotted with anti-MAST2 monoclonal antibody TA808790 (1:100). Then the blotted membrane was stripped and reprobed with anti-HSP90 antibody as a loading control.